E. Mettler  
Pelikanstrasse 19  
Zurich

---

*Mettler*

---

Switzerland


Suisse

Schweiz

---

Precision Balances  
Balances de précision  
Präzisionswaagen

---

K7/T  
K5/T  
K4/T

---

**Directions for use  
Mode d'emploi  
Betriebsvorschriften**

The Mettler Precision Balance is an instrument of practical and rugged design. It is simple to operate and in principle requires the observance of the following

### 3 Fundamental Rules for Maintaining Accuracy and Reliability

- 1 Set up balance where it will be free from vibration. A weighing table of adequate weight and stability will be a suitable base.
- 2 Lay the object to be weighed carefully onto the pan. Thus the knife-edges and bearing planes are subjected to less wear and tear.
- 3 The balance need no longer be arrested, between the individual weighing operations. When changing the location of the balance considerably, the weights should, however, be removed and the balance arrested.

english

K7/T

K5/T

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

## Unpacking

Unpack the balance and all the accessories carefully and clean them.

## Setting up


Fig. 1


Fig. 2

- 1 Release balance by turning knob ① (opposite to direction of arrow) up to stop. Check from time to time when using balance later on.
- 2 Remove cover ②.
- 3 Withdraw magnet protection ③.
- 4 Check voltage selector ④. The voltage shown must correspond to that of the mains. If necessary, readjust voltage (see «Adjustments», point 7.)
- 5 Turn knob ⑤ counter-clockwise till the figure 7 appears in the indicating window.
- 6 Clean the weights ⑥ and place them carefully into the forks ⑦.
- 7 Replace cover ②.

- 8 Slip pan carrier 8 onto the two guide pins 9.
- 9 Place pan 24 into position.
- 10 Insert power plug 10.

## Weighing


Fig. 1


Fig. 2

### 1 Checking the level

Check the spirit level 11 to ensure that balance is properly levelled. If necessary, relevel balance (see «Adjustments», point 1).

### 2 Checking zero point

- a) Unload pan and clean it if necessary.
- b) Turn knob 5 clockwise till the figure 0 appears in the indicating window.
- c) Switch on light with tumbler switch 12.
- d) **Only for K7/T:** Turn knob 13 (tare compensating device) clockwise up to stop.

K5/T

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

- e) With projected scale at rest, bring 0-line of scale to coincide with index line of screen by means of knob (14).
- f) If adjustment with knob (14) is no longer possible, proceed as under «Adjustments», point 2.

### 3 Laying on object to be weighed


Always place object to be weighed carefully on pan.

### 4 Weighing operations above 100 g

For weighing operations above 100 g, turn knob (5) counter-clockwise till the optical scale appears. Depending on the weight of the object to be weighed, the figures 1-7 will appear in the indicating window.

### 5 Indication of weight

Read the balance with the help of the index line (Fig. 3).

Weight = 237.78 g

Fig. 3

### 6 Weighing operations with a tare

By turning knob (13) (tare compensating device) tare weights up to 100 g can be balanced, i. e. the optical scale is reset to 0.

### 7 Weighing-in

By weighing-in operations, the approximate filling weight can be currently read from the point of intersection of the vertical oblique line with the small horizontal scale (below in the scale window).

## 8 Suspension lead-through

1. For balances which are supplied with suspension lead-through for weighings below balance.

With this type the object to be weighed can be secured to the hook below the base without special provisions.


Fig. 4

2. Rebuilding balances already supplied into type with suspension lead-through.
  - a) Arrest balance by turning knob ① in direction of arrow up to stop.
  - b) Free opening in base for the suspension lead-through by turning cover out of the way.
  - c) Remove hook ②5 from pan carrier ⑧, insert it through opening in base and hook it to suspension.
  - d) Release balance by turning knob ① (opposite to direction of arrow) up to stop.


Fig. 5


Fig. 6

K5/T

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

# Adjustments


Fig. 1


Fig. 2

Every balance is adjusted before leaving our works so that only minor readjustments may be needed as a result of shipping.

## 1 Levelling balance

Before doing any work on the balance, check whether it is properly levelled. If necessary, level the balance as follows:

- Observe spirit level (11) from above and adjust it in the longitudinal axis by means of the two front foot screws (15).
- Using the rear foot screw (15), bring spirit level (11) exactly into circle.

## 2 Zero point

If adjustment of the zero point is no longer possible with knob (14), proceed as follows:

**Only for K7/T:** Turn knob (13) (tare compensating device) clockwise up to stop.

- Level balance accurately (if balance was badly out of level, try once more to adjust zero point with knob (14)).
- Remove pan (24), pan carrier (8) and cover (2).

- c) Replace pan carrier (8) and pan (24).
- d) Adjust horizontal fine taring weight (16). Turning fine taring weight to the rear raises the 0-line. Turning fine taring weight to the front lowers the 0-line (previous to this, turn knob (14) back two full turns from end position).
- e) Remove pan (24) and pan carrier (8) replace cover (2) and then replace pan carrier (8) and pan (24).

**Caution:** Every zero point correction by means of the fine taring weight (16) must be followed by a check of the sensitivity (see points 3 a-h).

### 3 Sensitivity

The sensitivity of the balance can be checked and adjusted as follows:

- a) Remove pan (24), pan carrier (8) and cover (2).
- b) Replace pan carrier (8) and pan (24).
- c) Load pan with 100 g ( $\pm 0.2$  g). Turn knob (5) till the figure 1 appears in the indicating window.  
**Only for K7/T:** Turn knob (13) (tare compensating device) clockwise up to stop.
- d) Set zero point accurately.
- e) Carefully turn knob (5) clockwise till the figure 0 appears in indicating window.
- f) If the resultant deflection is smaller than 1000 graduations, turn adjusting weight (26) counter-clockwise and repeat operations c-e till a deflection of 1000 graduations is obtained.
- g) If the deflection is greater than 1000 graduations, turn adjusting weight (26) clockwise and repeat operations c-e till a deflection of 1000 graduations is obtained.
- h) Remove pan (24) and pan carrier (8), replace cover (2) and then replace pan carrier (8) and pan (24).

### 4 Adjusting definition of image

- a) A blurred scale image can be properly focused by turning the red screw (18).
- b) If quality of image deteriorates after long service, clean the mirrors (22) and (21) with enclosed brush.

### 5 Adjusting brilliance of image

- a) If the scale image is dark or shows coloured margins, readjust lamp (19) with knurled nut (20) provided on lampholder.
- b) See also point 4 b).

### 6 Displacing scale image laterally

The graduations of the moving scale should overlap the index line on the screen as shown in Fig. 3. If the overlap is too large or too small, light pressure on the mirror (21) in the direction of the arrow as per Fig. 1 will give the correct position.

K5/T

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

## 7 Setting to mains voltage

- a) Pull out power plug (10).
- b) Remove pan (24), pan carrier (8) and cover (2).
- c) Remove lamp (19).
- d) Pull out voltage selector knob (4). After turning it to required voltage, push it in again.
- e) Replace lamp (19).
- f) Replace cover (2), pan carrier (8) and pan (24).
- g) Plug in power plug (10).

## Maintenance


Fig. 1


Fig. 2

## Replacing lamp

After long service, it may happen that the lamp (19) fails or that a dark layer has formed on the inside of the bulb.

To change lamp (19), lightly depress bulb and turn it clockwise at the same time till the two short bayonet pins are released in the socket. In this position, withdraw the lamp. Fitting the spare lamp (23) kept in the cover of the housing is done in reverse sequence. Adjust brilliance of image as per «Adjustments», point 5.

**Important:** Procure another spare lamp immediately.

The Mettler Precision Balance is an instrument of practical and rugged design. It is simple to operate and, in principle, requires the observance of the following

### 3 Fundamental Rules for maintaining Accuracy and Reliability

- 1 Set up balance where it will be free from vibration. A weighing table of adequate weight and stability will be a suitable base.
- 2 Lay the object to be weighed carefully onto the pan. Thus the knife-edges and bearing planes will be subjected to less wear and tear.
- 3 The balance need no longer be arrested, between the individual weighing operations. When changing the location of the balance considerably, the weight should, however, be removed and the balance arrested.

K5/T

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

## Unpacking

Unpack the balance and all the accessories carefully and clean them.

## Setting up


Fig. 1


Fig. 2

- 1 Release balance by turning knob ① (opposite to direction of arrow) up to stop. Check from time to time when using balance later on.
- 2 Remove cover ②.
- 3 Withdraw magnet protection ③.
- 4 Check voltage selector ④. The voltage shown must correspond to that of the mains. If necessary, readjust voltage (see «Adjustments», point 7).
- 5 Turn lever ⑤ to horizontal position till the figure 0 appears in the indicating window.
- 6 Clean the weight ⑥ and place it carefully into fork ⑦.
- 7 Replace cover ②.

- 8 Slip pan carrier 8 onto the two guide pins 9.
- 9 Place pan 24 into position.
- 10 Insert power plug.

## Weighing


Fig. 1


Fig. 2

### 1 Checking the level

Check the spirit level 11 to ensure that balance is properly levelled. If necessary, relevel balance (see «Adjustments», point 1).

### 2 Checking zero point

- a) Unload pan and clean it if necessary.
- b) Turn lever 5 to horizontal position till the figure 0 appears in the indicating window.
- c) Switch on light with tumbler switch 12.
- d) **Only for K5/T:** Turn knob 13 (tare compensating device) clockwise up to stop.

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

- e) With projected scale at rest, bring 0-line of scale to coincide with index line of screen by means of knob (14).
- f) If adjustment with knob (14) is no longer possible, proceed as under «Adjustments», point 2.

### 3 Laying on object to be weighed


Always place object to be weighed carefully onto pan.

### 4 Weighing operations above 1000 g

For weighing operations above 1000 g, turn lever (5) to vertical position. The figure 1 appears in the indicating window.

### 5 Indication of weight

Read balance with the help of the index line (Fig. 3).

Weight = 1237.8 g

Fig. 3

### 6 Weighing operations with a tare

By turning knob (13) (tare compensating device), tare weights up to 500 g can be balanced, i. e. the optical scale is reset to 0.

### 7 Weighing-in

By weighing-in operations, the approximate filling weight can be currently read from the point of intersection of the vertical oblique line with the small horizontal scale (below in the scale window).

## 8 Suspension lead-through

1. For balances which are supplied with suspension lead-through for weighings below balance.

With this type, the object to be weighed can be secured to the hook below the base without special provisions.

Fig. 4

2. Rebuilding balances already supplied into type having suspension lead-through.

- Arrest balance by turning knob ① in direction of arrow up to stop.
- Free opening in base for the suspension lead-through by turning cover out of the way.
- Remove hook ②⑤ from pan carrier ⑧, insert it through opening in base and hook it to suspension.
- Release balance by turning knob ① (opposite to direction of arrow) up to stop.


Fig. 5


Fig. 6

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

## Adjustments


Fig. 1


Fig. 2

Every balance is adjusted before leaving our works so that only minor readjustments may be needed as a result of shipping.

### 1 Levelling balance

Before doing any work on the balance, check whether it is properly levelled. If necessary, level the balance as follows:

- Observe spirit level (11) from above and adjust it in the longitudinal axis by means of the two front foot screws (15).
- Using the rear foot screw (15), bring spirit level (11) exactly into circle.

### 2 Zero point

If adjustment of the zero point is no longer possible with knob (14), proceed as follows:

**Only for K5/T:** Turn knob (13) (tare compensating device) clockwise up to stop.

- Level balance accurately (if balance was badly out of level, try once more to adjust zero point with knob (14)).
- Remove pan (24), pan carrier (8) and cover (2).
- Replace pan carrier (8) and pan (24).

- d) Adjust fine taring weight (16). Turning fine taring weight to the rear raises the 0-line. Turning fine taring weight to the front lowers the 0-line (previous to this, turn knob (14) back two full turns from end position).
- e) Remove pan (24) and pan carrier (8), replace cover (2) and then replace pan carrier (8) and pan (24).

**Caution:** Every zero point correction by means of the fine taring weight (16) must be followed by a check of the sensitivity (see points 3 a-h).

### 3 Sensitivity

The sensitivity of the balance can be checked and adjusted as follows:

- a) Remove pan (24), pan carrier (8) and cover (2).
- b) Replace pan carrier (8) and pan (24).
- c) Load pan (24) with 1000 g ( $\pm 2$  g) and turn lever (5) to vertical position till the figure 1 appears in the indicating window.

**Only for K5/T:** Turn knob (13) (tare compensating device) clockwise up to stop.

- d) Set zero point accurately.
- e) Carefully turn lever (5) into horizontal position till the figure 0 appears in the indicating window.
- f) If the resultant deflection is smaller than 1000 graduations, turn spindle (17) clockwise and repeat operations c-e till a deflection of 1000 graduations is obtained.
- g) If deflection is greater than 1000 graduations, turn spindle (17) counter-clockwise and repeat operations c-e till a deflection of 1000 graduations is obtained.
- h) Remove pan (24) and pan carrier (8), replace cover (2) and then replace pan carrier (8) and pan (24).

### 4 Adjusting definition of image

- a) A blurred scale image can be properly focused by turning the red screw (18).
- b) If quality of image deteriorates after long service, clean the mirrors (22) and (21) with enclosed brush.

### 5 Adjusting brilliance of image

- a) If the scale image is dark or shows coloured margins, readjust lamp (19) with knurled nut (20) provided on lampholder.
- b) See also point 4 b).

### 6 Displacing scale image laterally

The graduations of the moving scale should overlap the index line on the screen as shown in Fig. 3. If the overlap is too large or too small, light pressure on the mirror (21) in the direction of the arrow as per Fig. 1 will give the correct position.

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

## 7 Setting to mains voltage

- a) Pull out power plug (10).
- b) Remove pan (24), pan carrier (8) and cover (2).
- c) Remove lamp (19).
- d) Pull out voltage selector knob (4). After turning it to required voltage, push it in again.
- e) Replace lamp (19).
- f) Replace cover (2), pan carrier (8) and pan (24).
- g) Plug in power plug (10).

## Maintenance


Fig. 1

Fig. 2

## Replacing lamp

After long service, it may happen that the lamp (19) fails or that a dark layer has formed on the inside of the bulb.

To change lamp (19), lightly depress bulb and turn it clockwise at the same time till the two short bayonet pins are released in the socket. In this position, withdraw the lamp. Fitting the spare lamp (23) kept in the cover of the housing is done in reverse sequence. Adjust brilliance of image as per «Adjustments», point 5.

**Important:** Procure another spare lamp immediately.

The Mettler Precision Balance is an instrument of practical and rugged design. It is simple to operate and, in principle, requires the observance of the following

### 3 Fundamentals Rules for Maintaining Accuracy and Reliability

- 1 Set up balance where it will be free from vibration. A weighing table of adequate weight and stability will be a suitable base.
- 2 Lay the object to be weighed carefully onto the pan. Thus the knife-edges and bearings planes will be subjected to less wear and tear.
- 3 The balance need no longer be arrested, between the individual weighing operations. When changing the location of the balance considerably, the weights should, however, be removed and the balance arrested (according to instructions as per "Weighing", point 8, 2a-2d).

K4/T

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

# Unpacking

Unpack the balance and all the accessories carefully and clean them.

## Setting up


Fig. 1


Fig. 2

- 1** Release balance by turning knob **1** (opposite to direction of arrow) up to stop. Check from time to time when using balance later on.
- 2** Remove cover **2**.
- 3** Withdraw magnet protection **3**.
- 4** Check voltage selector **4**. The voltage shown must correspond to that of the mains. If necessary, readjust voltage (see «Adjustments», point 7).
- 5** Turn knob **5** clockwise till the figure 0 appears in the indicating window.
- 6** Clean weight **6** and place it carefully into fork **7** with the help of pliers **27** (insert at an angle).

- 7 Clean weight (28) and carefully place it on flange (29).
- 8 Replace cover (2).
- 9 Slip pan carrier (8) onto the two guiding pins (9).
- 10 Place pan (24) into position.
- 11 Insert power plug.

## Weighing


Fig. 1


Fig. 2

### 1 Checking the level

Check the spirit level (11) to ensure that balance is properly levelled. If necessary, relevel balance (see «Adjustments», point 1).

### 2 Checking zero point

- a) Unload pan and clean it if necessary.
- b) Turn knob (5) clockwise till the figure 0 appears in the indicating window.
- c) Switch on light with tumbler switch (12).

- d) **Only for K4/T:** Turn knob ⑬ (tare compensating device) clockwise up to stop.
- e) With the projected scale at rest, bring 0-line of scale to coincide with index line of screen by means of knob ⑭.
- f) If adjustment with knob ⑭ is no longer possible, proceed as under «Adjustments», point 2.

### ③ Laying on object to be weighed


Always place object to be weighed carefully onto pan.

### ④ Weighing operations above 1000 g

For weighing operations above 1000 g, turn knob ⑤ counter-clockwise till optical scale appears. Depending on the weight of the object to be weighed, the figures 1-3 will appear.

### ⑤ Indication of weight

Read balance with the help of the index line (Fig. 3).

Weight = 3257.7 g

Fig. 3

### ⑥ Weighing operations with a tare

By turning the knob ⑬ (tare compensating device), tare weights up to 500 g can be balanced, i. e. the optical scale is reset to 0.

### ⑦ Weighing-in

By weighing-in operations, the approximate filling weight can be currently read from the point of intersection of the vertical oblique line with the small horizontal scale (below in the scale window).

## 8 Suspension lead-through

1. For balances which are supplied with suspension lead-through for weighings below balance.

With this type, the object to be weighed can be secured to the hook below the base without special provisions.

Fig. 4

2. Rebuilding balances already supplied into type having suspension lead-through.

- a) Remove pan (24), pan carrier (8) and cover (2).
- b) Turn knob (5) till the figure 3 appears in the indicating window.
- c) Remove weights (28) and (6).
- d) Arrest balance by turning knob (1) in direction of arrow up to stop.
- e) Free opening in base for the suspension lead-through by turning cover out of the way.
- f) Remove hook (25) from pan carrier (8), insert it through opening in base and hook it to suspension.
- g) Release balance in reversed sequence of operations as per points d-a.


Fig. 5


Fig. 6

français

K7/T

K5/T

K4/T

deutsch

K7/T

K5/T

K4/T

## Adjustments


Fig. 1


Fig. 2

Every balance is adjusted before leaving our works so that only minor readjustments may be needed as a result of shipping.

### 1 Levelling balance

Before doing any work on the balance, make sure it is properly levelled. If necessary, level the balance as follows:

- Observe spirit level (11) from above and adjust it in the longitudinal axis by means of the two front foot screws (15).
- Using the rear foot screw (12), bring the spirit level (11) exactly into circle.

### 2 Zero point

If adjustment of the zero point is no longer possible with the knob (14), proceed as follows:

**Only for K4/T:** Turn knob (13) (tare compensating device) clockwise up to stop.

- Level balance accurately (if balance was badly out of level try once more to adjust zero point with knob (14)).
- Remove pan (24), pan carrier (8) and cover (2).
- Replace pan carrier (8) and pan (24).

- d) Adjust fine taring weight (16). Turning fine taring weight to the rear raises the 0-line. Turning fine taring weight to the front lowers the 0-line (previous to this, turn knob (14) back two full turns from end position).
- e) Remove pan (24) and pan carrier (8), replace cover (2) and then replace pan carrier (8) and pan (24).

**Caution:** Every zero point correction by means of the fine taring weight (16) must be followed by a check of the sensitivity (see points 3 a-h).

### 3 Sensitivity

The sensitivity of the balance can be checked and adjusted as follows:

- a) Remove pan (24), pan carrier (8) and cover (2).
- b) Replace pan carrier (8) and pan (24).
- c) Load pan (24) with 1000 g ( $\pm 2$  g). Turn knob (5) till the figure 1 appears in the indicating window.  
**Only for K4/T:** Turn knob (13) (tare compensating device) clockwise up to stop.
- d) Set zero point accurately.
- e) Carefully turn knob (5) clockwise till the figure 0 appears in the indicating window.
- f) If the resultant deflection is smaller than 1000 graduations, turn spindle (17) clockwise and repeat operations c-e till a deflection of 1000 graduations is obtained.
- g) If the deflection is greater than 1000 graduations, turn spindle (17) counter-clockwise and repeat operations c-e till a deflection of 1000 graduations is obtained.
- h) Remove pan (24) and pan carrier (8), refit cover (2) and then replace pan carrier (8) and pan (24).

### 4 Adjusting definition of image

- a) A blurred scale image can be properly focused by turning the red screw (18).
- b) If quality of image deteriorates after long service, clean the mirrors (22) and (21) with enclosed brush.

### 5 Adjusting brilliance of image

- a) If the scale image is dark or shows coloured margins, readjust lamp (19) with knurled nut (20) provided on lamp holder.
- b) See also point 4 b).

### 6 Displacing scale image laterally


The graduations of the moving scale should overlap the index line on the screen as shown in Fig. 3. If the overlap is too large or too small, light pressure on the mirror (21) in the direction of the arrow as per Fig. 1 will give the correct position.

## 7 Setting to mains voltage

- a) Pull out power plug (10).
- b) Remove pan (24), pan carrier (8) and cover (2).
- c) Remove lamp (19).
- d) Pull out voltage selector knob (4). After turning it to required voltage, push it in again.
- e) Replace lamp (19).
- f) Replace cover (2), pan carrier (8) and pan (24).
- g) Plug in power plug (10).

## Maintenance

## Replacing lamp

After long service, it may happen that the lamp (19) fails or that a dark layer has formed on the inside of the bulb.

To change lamp (19), lightly depress bulb and turn it clockwise at the same time till the two short bayonet pins are released in the socket. In this position, withdraw the lamp. Fitting the spare lamp (23) kept in the cover of the housing is done in reverse sequence. Adjust brilliance of image as per «Adjustments», point 5.

**Important:** Procure another spare lamp immediately.

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